

IMPLANT CATALOG





IMPLANT CATALOG


Parlay Medical

Parlay Medical has supplied dental implants and dental implant products on a world-wide basis for many years.


Our product portfolio has four implant lines accompained by the richest set of accessories and tools in the industry with over 600 SKUs.

Our products are manufactured under quality management system that complies with the requirements of ISO 13485:2016 and are FDA approved, and are CE and MD Sap certified.





Five Internal Hex Implant Lengths





(Above renderings not to scale) Implants continued on next page...




Parlay Internal Hex Implants

are compatible with a number of internal hex implant systems such as Zimmer, Implant Direct, MIS, Hi-Tec and Bio Horizons














Internal Hex Implant Packaging




Abutment Options Shown Below








The Internal Hex Implant from Parlay Medical With 45 degree internal bevel






Wide-Pitch Thread Internal Hex Implants



Five Internal Hex Implant Lengths



(Above renderings not to scale)



The Single Platform (SP) Implant System

Five implant diameters with a 3.5mm platform:









3.2 3.7 4.2 4.7


5.2mm

- Platform switching for the 4.2, 4.7 and 5.2mm implants.
- Available with Narrow-Pitch (NP) or Wide-Pitch (WP) threads.

 Fixture mount is for implant insertion and can be used as a transfer and final straight abutment.


Standard Abutment






3.5mm Platform Color: Green





Wide-Pitch Implant



4.7 / 5.2 Implant Product Codes with 3.5 Platform:

NP: IHPXXXX


WP: IHWPXXXX

(Above renderings not to scale)



Internal Hex Implant Product Codes

		NP	THREAD	WP TH	READ
Diameter	Length	Cervical Contoured Abutment	Standard Abutment	Cervical Contoured Abutment	Standard Abutment
3.2mm (3.5	Platform)				
	8mm	IH328A	IH328	IHW328A	IHW328
	10mm	IH3210A	IH3210	IHW3210A	IHW3210
	11.5mm	IH3211A	IH3211	IHW3211A	IHW3211
	13mm	IH3213A	IH3213	IHW3213A	IHW3213
	16mm	IH3216A	IH3216	IHW3216A	IHW3216
3.7mm (3.5	Platform)				
	8mm	IH378A	IH378	IHW378A	IHW378
	10mm	IH3710A	IH3710	IHW3710A	IHW3710
	11.5mm	IH3711A	IH3711	IHW3711A	IHW3711
	13mm	IH3713A	IH3713	IHW3713A	IHW3713
	16mm	IH3716A	IH3716	IHW3716A	IHW3716
4.2mm (3.5	Platform)				
•	8mm	IH428A	IH428	IHW428A	IHW428
	10mm	IH4210A	IH4210	IHW4210A	IHW4210
	11.5mm	IH4211A	IH4211	IHW4211A	IHW4211
	13mm	IH4213A	IH4213	IHW4213A	IHW4213
	16mm	IH4216A	IH4216	IHW4216A	IHW4216
4.7mm (4.5	Platform)				
•	8mm	IH478A	IH478	IHW478A	IHW478
	10mm	IH4710A	IH4710	IHW4710A	IHW4710
	11.5mm	IH4711A	IH4711	IHW4711A	IHW4711
	13mm	IH4713A	IH4713	IHW4713A	IHW4713
	16mm	IH4716A	IH4716	IHW4716A	IHW4716
5.2mm (4.5	Platform)				
•	8mm	IH528A	IH528	IHW528A	IHW528
	10mm	IH5210A	IH5210	IHW5210A	IHW5210
	11.5mm	IH5211A	IH5211	IHW5211A	IHW5211
	13mm	IH5213A	IH5213	IHW5213A	IHW5213
	16mm	IH5216A	IH5216	IHW5216A	IHW5216
5.7mm (5.7	Platform)				
	8mm		IH578		
	10mm		IH5710		
	11.5mm		IH5711		
	13mm		IH5713		
	16mm		IH5716		



5 Internal Hex Implant Lengths

16mm, 13mm, 11.5mm, 10mm and 8mm

4.7 / 5.2 Implants (3.5 Platform)

NP	WP
IHP478	IHWP478
IHP4710	IHWP4710
IHP4711	IHWP4711
IHP4713	IHWP4713
IHP4716	IHWP4716
IHP528 IHP5210 IHP5211 IHP5213 IHP5216	IHWP5210 IHWP5211 IHWP5213 IHWP5216

NP + WP - Standard Mount



SURGICAL DRILLS

The step drills are coated to improve the visualization of the score lines and to increase cutting lubricity. They are available in 7 diameters. The first 4mm is the smaller diameter stepping up to the larger diameter.



Drill Sequence Chart								
			Platform/Im	plant Diamete	er			
Drill Code*	Size	Color	3.5/3.2	3.5/3.7	3.5/4.2	4.5/4.7	4.5/5.2	5.7/5.7
SD2320	2.3/2.0	White	X	X	X	X	X	X
SD2823	2.8/2.3	Grey	X	X	X	X	X	X
SD3428	3.4/2.8	Green		X	X	X	X	X
SD3834	3.8/3.4	Pink			X	X	X	X
SD4438	4.4/3.8	Purple				X	X	X
SD4844	4.8/4.4	Blue					X	X
SD5448	5.4/4.8	Yellow						X


Last drill is used in dense bone situations. For soft bone, delete last drill.

(Above photos not to scale)

Implant Instrumentation continued on next page...

^{*} For use with drill stops.





Precise Adjustable Torque Wrench with Nine Settings: 15, 20, 25, 30, 35, 40, 45, 50 and 60 Ncm

ATW




Single Setting Torque Wrenches: 15, 20, 25, 30 and 35 Ncm 15Ncm: GTW15; 20Ncm: GTW20; 25Ncm: GTW25; 30Ncm: GTW30; 35Ncm: GTW35


(Above photos not to scale)

Implant Instrumentation continued on next page...


10









4mm Square Screw Driver Handle: ISIH




Implant Mount Driver-Short: IMDSN Long: IMDLN



4mm Ratchet Wrench: RWA



Drill Extension: DEXI



Implant Instrumentation continued on next page...

(Above photos not to scale)



Implant Drivers, Extensions & Adaptors

SHAFT LENGTHS: SHORT: 6.4mm

SPECIALTY DRIVERS SPECS .048 /.050 Torx NobelBiocare Torx Short: ITIS Long: ITIL 4mm Square Driver Short: TRIS 4mm Square Driver: Long: TRIL .080 Hex Short: HFD **TGDN** Mini Implant Driver Short MIDS Mini Implant




(Product photos not to scale)




Implant Drivers, Extensions & Adaptors







Bone Expander Kit



An Alternative to Osteotomes

Bone expander drills are an alternative to osteotomes for the expansion and condensing of the atrophic mandible and maxilla in preparation for dental implant insertions. Expanders are also an alternative to the maxillary sinus elevation technique.

Surgical Drama Decreased

Expanders are driven into the bone with a ratchet wrench or low speed hand piece. This decreases the surgical trauma of osteotomes. Bone expanders improve the clinical success by improving stability, maintaining bone density and increasing fixation.

What the Kit Includes

The Bone Expander Kit from Parlay Medical includes either 5 or 6 bone expanders, a ratchet wrench, a pilot drill, a ratchet wrench extender, a hand piece latch adaptor and a manual knob for finger usage. The kit is autoclavable.

Product code: BEK

14



Basic Surgical Drilling Kit



The Parlay Medical Basic Surgical Drilling Kit Product Code: SKK

The Surgical Drilling Kit shown here includes a Ratchet Wrench, 1.25mm Hex Driver, Motor Driver, Manual Knob,

Pilot Drill, 2.5mm Hex Driver, Long, 2.5mm Hex Driver, Short, 2 Paralleling Pins, with the following Drills: 2.3/2.0; 2.8/2.3; 3.4/2.8; 3.8/3.4; 4.4/3.8; 4.8/4.4 and 5.4/4.8mm Drills. The Kit is fully autoclaveable and constructed of high-grade stainless steel.

Autoclavable:

This kit may be autoclaved. Follow the sterilizer manufacturer's recommended procedure while observing the following cautions:

Caution:

- Do not exceed 275° (175°C)
- Do not submerge in any solution
- Do not use ultrasonic devices
- **Do not** autoclave for extended periods of time (i.e. overnight). Extended autoclaving should not be confused with repeated autoclaving. Repeated autoclaving will not harm the product.

(Above photos not to scale)

Implant Instrumentation continued on next page...





Advanced Surgical Drilling Kit



The Advanced Surgical

Drilling Kit shown here includes:

Ratchet Wrench, Straight Handle, Drill Extension, Implant Mount Driver, Thumb Knob, Long and Short .048 Hex Driver, Long and Short .050 Hex Driver, Long and Short 2.5mm Hex Driver, 4, Direction Indicators, Pilot Drill, and 7, Surgical Drills: 2.3/2.0, 2.8/2.3, 3.4/2.8, 3.8/3.4, 4.4/3.8, 4.8/4.4, and 5.4/4.8. The Kit is autoclaveable.

Autoclaving:

This kit may be autoclaved. Follow the sterilizer manufacturer's recommended procedure while observing the following cautions:

Caution:

- Do not exceed 275° (175°C)
- Do not submerge in any solution
- Do not use ultrasonic devices
- Do not autoclave for extended periods of time (i.e. overnight). Extended autoclaving should not be confused with repeated autoclaving. Repeated autoclaving will not harm the product.




Prosthetic Catalog Index Internal Hex

Implant Items	IMPLANT ITEMS AND INSTRUMENTATION	I - 15
Section I	OLS SYSTEM	18-19
Section 2	MLS SYSTEM	20-21
Section 3	MLS UNIT RESTORATION	22-23
Section 4	BALL ATTACHMENTS	24-26
Section 5	TEMERORARY ABUTMENTS	27
Section 6	TITANIUM ABUTMENTS	28-3 I
Section 7	CASTING ABUTMENTS	32
Section 8	ANALOG	33
Section 9	3D SCANNING INTERFACE	34
Section 10	HEALING CAPS	35
Section I I	PROSTHETIC SCREWS	36
Section 12	TRANSFER	37



Prosthetic Catalog OLS System





OLS ABUTMENT SYSTEM - ACCESSORIES






Prosthetic Catalog
OLS System


ANGULAR OLS ABUTMENT SYSTEM CONCEPT 9°,18°,30°


OLS Over Denture Attachment - SETS


9° OLS ABUTMENT SYSTEM ANGULAR MLS SYSTEM CONCEPT

OLS Angular Over Denture Attachment - SET










Cat. No.

Designation

Cat. No.

Designation



30° OLS ABUTMENT SYSTEM ANGULAR MLS SYSTEM CONCEPT


MLS System 3 mm

Height

OLS Angular Over Denture Attachment - SET

MLS System 2 mm

Height



Designation info for Angular MLS System Adaptors:

Internal Hex:

3.5 platform = WCA-... - M

4.5 platform = WCA-... - W

MLS System 5 mm

Height

MLS System 4 mm


Height



Prosthetic Catalog MLS System

ANGULAR MLS SYSTEM CONCEPT 9°,18°,30° MLS Unit Items Over Denture Attachment

MLS SYSTEM





9° ANGLE MLS

SYSTEM Adaptor

Cat. No.

Designation



18° ANGLE MLS

SYSTEM Adaptor



30° ANGLE MLS

SYSTEM Adaptor



WIDE UNIT - System Adaptor Angular MLS System Adaptor 9°,18°,30° 4.5 Platform







Cat. No.	EWU-0910	EWU-1810	EWU-3010
Designation	9° ANGLE MLS	18° ANGLE MLS	30° ANGLE MLS
	SYSTEM Adaptor	SYSTEM Adaptor	SYSTEM Adaptor

ANGLE AND HEIGHT GAUGE












Cat. No.	GAHM-00	GAHM-09	GAHM-18	GAHM-30	GAHM-SET
Designation	Straight gauge	9° Angle gauge	18° Angle gauge	30° Angle gauge	Set include full range of gauges

SYSTEM CONCEPT





WCC-0020

2 mm CONNECTION

MLS SYSTEM

OLS



WCC-0030

3 mm CONNECTION

MLS SYSTEM

OLS



WCC-0040

4 mm CONNECTION

MLS SYSTEM

OLS



WCC-0050

5 mm CONNECTION

MLS SYSTEM

OLS






OLS ABUTMENT SYSTEM - ACCESSORIES

Cat. No.

Designation











WC-1301	
NALOC for	

Cat. No.	WCH-5628	WCH-4823	WCP-2348-H	WCPS-2348-H	WCP-2348	SCH-8540	TWC-1105	WCT-900	AWC-1301
Designation	STAINLESS STEEL HOUSING COVER	STAINLESS STEEL HOUSING	Hard Housing PEEK Attachment	Soft Housing PEEK Attachment	DIRECT PEEK Attachment	SILICON Sealed RING	TRANSFER for OLS	OLS Tool	ANALOG for OLS



Prosthetic Catalog MLS System

BALL MLS SYSTEM ANGULAR

Ball attachment connection for MLS System Concept



BMA-0020

2 mm CONNECTION

MLS SYSTEM BALL







BMA-0030	BMA-0040	BMA-0050
3 mm CONNECTION	4 mm CONNECTION	5 mm CONNECTION
MLS SYSTEM BALL	MLS SYSTEM BALL	MLS SYSTEM BALL

BALL Attachment- Accessories





Cat. No.	BHM-5032	BIS-0025	
Designation	Stainless Steel Housing	Silicon Insert	



Cat. No.

Designation

2 MLS SYSTEM COI

MLS SYSTEM UNIT CONNECTION

Screwed in Connection for MLS System Concept



MUC-0021

2 mm CONNECTION

MLS SYSTEM



MUC-0031

3 mm CONNECTION

MLS SYSTEM




MUC-0041

4 mm CONNECTION

MLS SYSTEM



MLS SYSTEM



CONCEPT

Color no. 2 Internal Hex Cat. No.

Designation

MLS-UNIT ACCESORIES





NOTE: It is recommended to purchase additional screws.



Prosthetic Catalog MLS Unit Restoration

MLS SYSTEM STRAIGHT AESTHETIC ABUTMENTS

Abutments with smooth, non locking surface for MLS-Unit restorations

AESTHETIC MINOR UNIT Straight MLS-Unit 3.5 platform





Designation












Cat. No.	EMU-0011	EMU-0021	EMU-0031	EMU-0041	EMU-0051	EMU-0061
Designation	1mm height	2mm height	3mm height	4mm height	5mm height	6mm height




Cat, No.

Designation



2mm height



MLS-UNIT ACCESORIES

1mm height

	PSU-0032 Plastic Sleeve screw included	PSU-0050 Parallel connection	TSU-1250 Titanium Sleeve screw included L = 12.0mm	TMU-0532 Titanium Short Sleeve srew included	





Titanium Sleeve

3mm height



Titanium Sleeve

L = 80mm

4mm height



3-D Scanning

Device



TMU-1000

Open Tray

Short Transfer

Cat. No.

Designation

5mm height



TMU-1200

Open Tray


Transfer

screw included









TMU-8005 Closed Tray Transfer

HCC-4550 WQP-0011 Plastic for Closed Healing Cap Tray Transfer for Multi-Unit

AMU-1201 Analog for Multi-Unit

OVER DENTURE ABUTMENTS

For single tooth restotations

Over Denture Abutments with Hex Locking for single tooth restorations.

AESTHETIC CONNECTION ABUTMENT with hex for Retention and positioning








Cat. No.	ECA-6005	ECA-6015	ECA-6025
Designation	0.5mm height	1.5mm height	2.5mm height
Related screw	SCS - 3011	SCS - 3012	SCS - 3013

PLASTIC AESTHETIC SLEEVE











Prosthetic Catalog MLS Unit Restoration






Prosthetic Catalog Ball Attachments

Over Denture Ball Attachment in Connection with Implants







BALL HOUSING MINOR / BALL INSERT SILICON


Cat. No.	BHM-5032	BIS-0025
Designation	Stainless housing	Silicon insert




ATTACHMENT

Color no. 4 Internal Hex

Prosthetic Catalog Ball Attachments





BALL HOUSING MINOR / BALL INSERT SILICON

Cat. No.	BHM-5032	BIS-0025
Designation	Stainless housing	Silicon insert




Prosthetic Catalog Ball Attachments

ANGULAR MLS SYSTEM CONCEPT 9°,18°,30°

Angular Ball Attachment - Set

BALL MINOR ANGULAR MLS SYSTEM CONCEPT


3.5 Platform



BALL MINOR ANGULAR MLS SYSTEM CONCEPT

3.5 Platform

Cat. No.




Designation info for **Angular MLS System Adaptors:**

Internal Hex: 3.5 Platform = BMA-...-S**4.5 Platform = BWA-....-S**

30° BALL MINOR ANGULAR MLS SYSTEM CONCEPT


3.5 Platform

ATTACHMENT



BALL HOUSING MINOR / BALL INSERT SILICON





OLH 0035 Cat. No. Orienting Lead Designation Holder



Prosthetic Catalog Temporary Abutments

TEMPORARY AESTHETIC ABUTMENTS

Includes a Retaining Screw

STRAIGHT LEVEL MINOR - PEEK 3.5 Platform







Cat. No.	SLM-1110P	SLM-1220P	SLM-1330P
Designation	1mm shoulder	2mm shoulder	2.5mm shoulder

Abutments for temporary restoration made of biocompatible material PEEK. Polyether ether Ether Ketone is a colorless organic Thermoplastic polymer, proven to provide long term service.


TEMPORARY TITANIUM ABUTMENT 3.5 Platform

Cat. No.

Designation

Cat. No.

Designation



Engaging





ANGULAR CURVE MINOR - PEEK Platform 15° 3.5



ACM-1115-10P

1mm Aesthetic

shoulder







ACM-1215-20P ACM-1215-30P ACM-1315-40P 2mm Aesthetic 4mm Aesthetic 3mm Aesthetic shoulder shoulder shoulder

IMMEDIATE TEMPORARY CONICAL



금

MPORARY ABUTMENTS


Color no.

Internal Hex












ACM-1225-20P 2mm Aesthetic shoulder








Prosthetic Catalog Titanium Abutments

STRAIGHT AESTHETIC ABUTMENTS

Includes a Retaining Screw

Ultra Slim Abutments are designed to be used in cases where the gaps between the implants are too narrow for standard abutments.







Prosthetic Catalog Titanium Abutments

ANGULAR AESTHETIC ABUTMENTS

Includes a Retaining Screw

SPECIFICATIONS

- The angular aesthetic abutment body functions as impression coping
- Includes the Prosthetic Screw
- Installs with 0.050 hex drivers (1.27mm)
- Non-rotational design for standard, internal-hex implants
- Made of Commercial grade 5 Titanium alloy (ASTM-F-136 02A COND A)



ANGULAR CURVE MINOR 3.5 Platform 15°









shoulder

0

ITANIUM ABUTMENTS

Cat. No. ACM-1115-10 ACM-1315-40 ACM-1215-20 ACM-1215-30 1mm Aesthetic 4mm Aesthetic 2mm Aesthetic 3mm Aesthetic Designation shoulder shoulder shoulder

ANGULAR CURVE WIDE

4.5 Platform 15°






Color

Internal Hex

Cat. No.	ACW-1115-10	ACW-1215-20
Designation	1mm Aesthetic shoulder	2mm Aesthetic shoulder



ANGULAR CURVE MINOR

3.5 Platform 25°

Cat. No.

Designation



ACM-1125-10

1mm Aesthetic


shoulder







ACM-1225-20	ACM-1325-30	ACM-1425-40
2mm Aesthetic shoulder	3mm Aesthetic sholuder	4mm Aesthetic shoulder





Prosthetic Catalog Cervical Contour Abutments

STRAIGHT & ANGLED CERVICAL CONTOURED ABUTMENTS

Includes a Retaining Screw

SPECIFICATIONS

- Installs with 0.050 hex drivers (1.27mm)
- Non-rotational design for standard, internal-hex implants
- Made of Commercial grade 5 Titanium alloy (ASTM-F-136 02A COND A)



STRAIGHT CERVICAL CONTOUR

3.5 Platform







Cat. No.	SCCM-1220
Designation	Straight Aesthetic Abutment 2mmH

SCCM-1330 Straight Aesthetic Abutment 3mmH SCCM-1440
Straight Aesthetic
Abutment 4mmH

6 CERVICAL CONTOUR ABUTMENTS

STRAIGHT CERVICAL CONTOUR

4.5 Platform







Cat. No.	SCCW-1220	SCCW-1330	SCCW-1440
Designation	Straight Aesthetic Abutment 2mmH	Straight Aesthetic Abutment 3mmH	Straight Aesthetic Abutment 4mmH

ANGULAR CERVICAL CONTOUR

3.5 Platform 15°

Cat. No.

Designation




15° Aesthetic

Abutment 2mmH







ANGULAR CERVICAL CONTOUR

3.5 Platform 25°







Cat. No.	ACCM-1225-20	ACCM-1325-30	ACCM-1425-40
Designation	25° Aesthetic Abutment 2mmH	25° Aesthetic Abutment 3mmH	25° Aesthetic Abutment 4mmH



Prosthetic Catalog Titanium Abutments

ANGULAR SMOOTH MINOR - LEAF

3.5 Platform 15°





Cat. No.	ASM-0915 - LF	ASM-1115 - LF
Designation	Chopped edge	Chopped edge

ANGULAR SMOOTH MINOR - LEAF 3.5 Platform 25°

Ultra Slim and LEAF treatment are for front teeth



ITANIUM ABUTMENTS

Color no.

Internal Hex





at. No.	ASM-0925 - LF	ASM-1125 - LF
esignation	Chopped edge	Chopped edge

ANGULAR SMOOTH MINOR 3.5 Platform 15°



ASM-0915

Smooth shape



ASM-1115

Smooth shape



W

ASM-1315
Smooth shape

ANGULAR SMOOTH THIN 3.5 Platform 15° & 25°



Cat. No.	AST-1115	AST-1125	
Designation	Ultra Slim 15°	Ultra Slim 25°	

ANGULAR SMOOTH WIDE 4.5 Platform 15°



Cat. No.	ASW-1115
Designation	Smooth shape

ANGULAR SMOOTH MINOR

Cat. No.

Designation

3.5 Platform 25°







Cat. No.	ASM-0925	ASM-1125	ASM-1325
Designation	Smooth shape	Smooth shape	Smooth shape




Prosthetic Catalog Casting Abutments

PLASTIC ABUTMENTS FOR CASTING

PLASTICS STRAIGHT FOR CASTING

3.5 / 4.5 Platform



ASTING **ABUTMENTS**

The Titanium / Plastic Abutments are for Investment and casting following conventional crown and bridge techniques.

The Titanium / Plastic Abutments combine precision machined interface between the implant and the Abutment with the convenience of a castable plastic sleeve.

TITANIUM PLASTIC ABUTMENT

3.5 Platform



Cat. No.

Designation





Cat. No.	TPA-1420 W	TPA-1420 WR
Designation	Hex Base wide platform	Round Base wide platform

CERVICAL CONTOURED PLASTIC ABUTMENTS

3.5 Platform, 3mm Shoulder






Cat. No.	CPA-1420	CPA-1420-R
Designation	Hex Base	Round Base





CPA-1430	CPA-1430-R
Hex Base H= 3.0 mm	Round Base H= 3.0 mm






Cat. No.	CPA-1420-W	CPA-1420-WR
Designation	Hex Base Wide Platform	Round Base Wide Platform




Prosthetic Catalog Analog

Analogs are used for the fabrication of a master model. They are replicas of the implant. All analogs are made of high grade stainless steel.

ANALOG STAINLESS










MLS SYSTEM

ANALOG FOR OLS ABUTMENT SYSTEM



OLS ABUTMENT SYSTEM



Prosthetic Catalog 3D Scanning Interface



3.5 Platform - Non-Rotational

SMS - 6031	SMS - 6032	SMS - 6033	SMS - 6034
H=01 CadCam	H=02 CadCam	H=03 CadCam	H=04 CadCam
TiBase	TiBase	TiBase	TiBase
Non-rotational	Non-rotational	Non-rotational	Non-rotational

3.5 Platform - Rotational



4.5 Platform

SWS - 6030	SWS - 0030		
FB			
CadCam TiBase wide	CadCam TiBase wide		
Non-rotational	Rotational		

NOTE: Use screw SCS - 3095 or SCP - 2580 for scanning abutments.



Prosthetic Catalog Healing Caps

Healing caps are used to maintain access through tissue during the healing phase. They are available in variety of heights and diameters to accommodate tissue depth.

HEALING CAP MINOR

3.5 Ultra Slim Platform













Cat. No.	HCM-3720	HCM-3730	HCM-3740	HCM-3750	HCM-3760	HCM-3770
Designation	Dia. 3.7 mm					
	2mm height	3mm height	4mm height	5mm height	6mm height	7mm height

^{*} Recomended for narrow spaces

HEALING CAP MINOR

3.5 Slim Platform



HEALING













Cat. No.	HCM-4020	HCM-4030	HCM-4040	HCM-4050	HCM-4060	HCM-4070
Designation	Dia. 4.6 mm					
	2mm height	3mm height	4mm height	5mm height	6mm height	7mm height

^{*} Recomended for narrow spaces

HEALING CAP MINOR 3.5 Normal Platform













Cat. No.	HCM-5020	HCM-5030	HCM-5040	HCM-5050	HCM-5060	HCM-5070
Designation	Dia. 5 mm 2mm height	Dia. 5 mm 3mm height	Dia. 5 mm 4mm height	Dia. 5 mm 5mm height	Dia. 5 mm 6mm height	Dia. 5 mm 7mm height

^{*} Recomended for angular restoration

HEALING CAP WIDE 4.5 Wide Platform

NOTE: Healing Caps with straight walls instead of the cervical contoured shown, are available and have an "S" at the end of the Product Code,

i.e.: HCM-5030S.














HCW-6080 HCW-6070 Dia. 6 mm Dia. 6 mm 8mm height 7mm height





Prosthetic Catalog Prosthetic Screws

SCREWS FOR ABUTMENTS



SCREWS FOR OVER DENTURE ABUTMENTS











Prosthetic Catalog Transfer

TRANSFER IMPRESSION COPING


Includes a Retaining Screw







* Plastic cap included



3.5 Platform



shoulder

shoulder

shoulder

shoulder

shoulder

Cat. No.	TMC-0900		
Designation			





Parlay Medical 617 N. Oxford Ave. #208 Los Angeles, CA 90004 +1-555-ParlayM (855-727-5296) www.parlaymedical.com







